



The Reliable and Rapid


DxTROP I

Troponin I card test

CHEST PAIN



SUSPECTED MI



**When Every
Minute Counts...!!**




**Point of care. Anytime.
Anywhere.**

DiaSys

Diagnostic Systems

CHOOSING QUALITY.



Highly Cardiac Specific Biomarker with cut - off level of 0.3 ng/mL

Troponin I is proven to be more reliable, sensitive & specific than Troponin T, Myoglobin & CKMB

- The use of a sensitive assay for Troponin I (Trop I) improves early diagnosis of acute myocardial infarction and risk stratification, regardless of the time of chest-pain onset.¹
- Trop I - specific for myocardial injury.²
- Cardiac Troponin I (cTnI) should be the biomarker of choice for diagnosis of myocardial injury.³
- cTnI has been shown to be 98.9% specific for the heart.⁴

Reference:
The New England journal of Medicine: Volume 361:868-877 August 27, 2009 Number 9


ACCET image provided by the American College of Cardiology Foundation.

Cardiac Enzymes, Renal Failure and Renal Transplantation Huseyin Bozbas, Aylin Yildirim, and Haldun Muderrisoglu
CLINICAL MEDICINE AND RESEARCH 2006 4:79-84

Department of Laboratory Medicine and Pathology, Hennepin County Medical Center, 701 Park Ave., Minneapolis, MN 55415, USA.

Testing Procedure and Interpretation of results

Step 1



Using dropper, transfer 1 drop (50µL) of whole blood or serum or plasma specimen into the sample well

Step 2




Add 1 drop of QDx Trop I Buffer

Step 3




Read the test result in 10 minutes




Positive result

If pink band appears in Control as well as in Test window, test is positive indicating high risk of myocardial infarction



Negative result

If pink appears only in control window then test negative



Invalid

If no pink band is present in the control window at 10 minutes, the test is invalid




Features

- Rapid & reliable one-step assay enabling fast results Test
- Results at just 10 minutes.
- Sensitivity: 96.9% & Specificity: 98.9%
- Whole blood/Serum/Plasma as sample. (Universal Sampling)
- Storage temperature at 2°-30°C



Kit Specifications	
Test Type	Rapid Qualitative test
Sample Type	Whole Blood/ Serum/ Plasma
Sample Size	1 drops
Detecting	Cardiac Troponin I (CTnI)
Test Principle	Lateral Flow Sandwich Immunoassay
Read time	10 minutes
Cut off value	0.3 ng/mL
Sensitivity	96.9%
Specificity	98.9%
Storage temperature	2° - 30°C

Kit Specifications	
1. QDx Trop I (10 tests)	Cat No : 400295508
Each box contain :	
a) 10 tests device with a desiccant pouch and a plastic dropper in each test.	
b) Pack insert : 1	
c) Assay buffer bottle : 1	
Pack insert : 1	
Assay buffer bottle : 1	
2. QDx Trop I (10 tests : Mono Carton)	Cat No : 400295548
Each mono carton contain :	
a) One test device with 1 desiccant pouch and 1 plastic dropper.	
b) Pack insert : 1	
c) Assay buffer bottle : 1	
d) Report	



DiaSys
Diagnostics India Private Limited
 Plot No. A -821,TTC, MIDC, Mahape,
 Navi-Mumbai, 400710.
 Maharashtra, India.

Toll free: 1800 120 1447
 E-Mail : info@diasys.in
www.diasys.in



Oct 2018

CHOOSING QUALITY.