

ASAT (GOT) FS* (IFCC mod.)

with/without pyridoxal-5-phosphate

Diagnostic reagent for quantitative in vitro determination of ASAT (GOT) in serum or plasma on DiaSys respons[®]910

Order Information

Cat. No. 1 2601 99 10 920

4 twin containers for 200 tests each

Pyridoxal-5-phosphate FS


Cat. No. 2 5010 99 10 030

6 x 3 mL

Method

Optimized UV-test according to IFCC (International Federation of Clinical Chemistry and Laboratory Medicine) [modified]

Principle



Addition of pyridoxal-5-phosphate (P-5-P) stabilizes the activity of transaminases and avoids falsely low values in samples containing insufficient endogenous P-5-P, e.g. from patients with myocardial infarction, liver disease and intensive care patients [1].

Reagents

Components and Concentrations

R1:	TRIS	pH 7.65	110 mmol/L
	L-Aspartate		320 mmol/L
	MDH (malate dehydrogenase)		≥ 800 U/L
	LDH (lactate dehydrogenase)		≥ 1200 U/L
R2:	2-Oxoglutarate		65 mmol/L
	NADH		1 mmol/L

Pyridoxal-5-phosphate FS

Good's buffer	pH 9.6	100 mmol/L
Pyridoxal-5-phosphate		13 mmol/L

Storage Instructions and Reagent Stability

The reagents are stable up to the end of the indicated month of expiry, if stored at 2 – 8 °C, protected from light and contamination is avoided. DiaSys respons containers provide protection from light. Do not freeze the reagents!

Warnings and Precautions

- The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
- In very rare cases, samples of patients with gammopathy might give falsified results.
- Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examinations and other findings

Waste Management

Please refer to local legal requirements.

Reagent Preparation

The reagents are ready to use. The bottles are placed directly into the reagent rotor. For the determination with P-5-P add 350 µL of P-5-P to reagent 1 and mix gently.

Stability after mixing: 6 days at 2 - 8 °C
24 hours at 15 - 25 °C

Specimen

Serum, heparin plasma or EDTA plasma

Stability [2]:

4 days at 20 - 25 °C

7 days at 4 - 8 °C

3 months at -20 °C

Discard contaminated specimens. Freeze only once.

Calibrators and Controls

For the calibration, DiaSys TruCal U calibrator is recommended. This method has been standardized against the original IFCC formulation. For internal quality control DiaSys TruLab N and P controls should be assayed. Each laboratory should establish corrective action in case of deviations in control recovery.

	Cat. No.	Kit size
TruCal U	5 9100 99 10 063	20 x 3 mL
	5 9100 99 10 064	6 x 3 mL
TruLab N	5 9000 99 10 062	20 x 5 mL
	5 9000 99 10 061	6 x 5 mL
TruLab P	5 9050 99 10 062	20 x 5 mL
	5 9050 99 10 061	6 x 5 mL

Performance Characteristics

with P-5-P activation

Measuring range up to 675 U/L ASAT (in case of higher activities re-measure samples after manual dilution or use rerun function)	
Limit of detection**	2 U/L ASAT
On-board stability	6 days
Calibration stability	6 days

Interfering substance	Interferences < 10%	ASAT [U/L]
Ascorbate	up to 30 mg/dL	108
Hemoglobin	up to 20 mg/dL	22.9
Bilirubin, conjugated	up to 100 mg/dL	166
	up to 55 mg/dL	42.6
Bilirubin, unconjugated	up to 55 mg/dL	165
	up to 60 mg/dL	44.0
	up to 60 mg/dL	173
Lipemia (triglycerides)	up to 1000 mg/dL	39.2
	up to 500 mg/dL	149

For further information on interfering substances refer to Young DS [4].

Precision			
Within run (n=20)	Sample 1	Sample 2	Sample 3
Mean [U/L]	35.1	44.4	172
Coefficient of variation [%]	1.54	1.85	1.47
Between run (n=20)	Sample 1	Sample 2	Sample 3
Mean [U/L]	27.9	44.7	174
Coefficient of variation [%]	4.07	2.71	1.34

Method comparison (n=115)	
Test x	DiaSys ASAT (GOT) FS (Hitachi 917)
Test y	DiaSys ASAT (GOT) FS (respons [®] 920)
Slope	1.033
Intercept	-2.309 U/L
Coefficient of correlation	0.9996

** Lowest measurable concentration which can be distinguished from zero mean + 3 SD (n = 20) of an analyte free specimen

without P-5-P activation

Measuring range up to 700 U/L ASAT (in case of higher activities re-measure samples after manual dilution or use rerun function)	
Limit of detection***	2 U/L ASAT
Onboard stability	4 weeks
Calibration stability	4 weeks

Interfering substance	Interferences < 10%	ASAT [U/L]
Ascorbate	up to 30 mg/dL	125
Hemoglobin	up to 50 mg/dL	22.6
Bilirubin, conjugated	up to 10 mg/dL	19.0
	up to 65 mg/dL	36.7
Bilirubin, unconjugated	up to 70 mg/dL	18.6
Lipemia (triglycerides)	up to 1000 mg/dL	43.7
	up to 1300 mg/dL	175
For further information on interfering substances refer to Young S [4].		

Precision			
Within run (n=20)	Sample 1	Sample 2	Sample 3
Mean [U/L]	23.5	40.1	199
Coefficient of variation [%]	2.54	1.61	1.07
Between run (n=20)	Sample 1	Sample 2	Sample 3
Mean [U/L]	25.5	49.4	205
Coefficient of variation [%]	3.13	1.55	1.00

Method comparison (n=105)	
Test x	DiaSys ASAT (GOT) FS (Hitachi 917)
Test y	DiaSys ASAT (GOT) FS (respons [®] 910)
Slope	0.996
Intercept	0.079 U/L
Coefficient of correlation	0.9999

*** Lowest measurable concentration which can be distinguished from zero mean + 3 SD (n = 20) of an analyte free specimen

Conversion factor

ASAT [U/L] x 0.0167 = ASAT [μkat/L]

Reference Range

With pyridoxal-5-phosphate activation

Women [3]	< 31 U/L	< 0,52 μkat/L
Men [3]	< 35 U/L	< 0,58 μkat/L
Children [1]	1 – 3 year(s)	< 50 U/L
	4 – 6 years	< 45 U/L
	7 – 9 years	< 40 U/L
	10 – 12 years	< 40 U/L
	13 – 15 years	< 35 U/L
	16 – 18 years	< 35 U/L
		< 0,83 μkat/L
		< 0,75 μkat/L
		< 0,67 μkat/L
		< 0,67 μkat/L
		< 0,58 μkat/L
		< 0,58 μkat/L

Without pyridoxal-5-phosphate activation

Women	< 31 U/L	< 0,52 μkat/L
Men	< 35 U/L	< 0,58 μkat/L

Each laboratory should check if the reference ranges are transferable to its own patient population and determine own reference ranges if necessary.

Literature

1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 18-9.
3. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 5: Reference procedure for the measurement of catalytic concentration of aspartate aminotransferase. Clin Chem Lab Med 2002; 40: 725-33.
4. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th. ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press, 2000.
5. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.

Manufacturer



DiaSys Diagnostic Systems GmbH
Alte Strasse 9 65558 Holzheim Germany

ASAT (GOT) FS (IFCC mod.)

Application for serum and plasma samples

This application was set up and evaluated by DiaSys. It is based on the standard equipment at that time and does not apply to any equipment modifications undertaken by unqualified personnel

Identification	
This method is usable for analysis:	Yes
Name:	AST
Shortcut:	
Reagent barcode reference:	011
Host reference:	

Technic	
Type:	Linear Kinetic
First reagent:[μ L]	160
Blanc correction	Yes
Second reagent:[μ L]	40
Blanc correction	Yes
Main wavelength:[nm]	340
Secondary wavelength:[nm]	405
Polychromatic factor:	1.000
1 st reading time [min:sec]	06:24
Last reading time [min:sec]	10:00
Reaction way:	Decreasing
Linear Kinetics	
Substrate depletion: absorbance limit	0.3
Linearity: Maximum deviation [%]	100
Fixed Time Kinetics	
Substrate depletion: absorbance limit	
Endpoint	
Stability: largest remaining slope	
Prozone Limit [%]	

Sample	
Diluent	NaCl
Concentration technical limits-Lower	2.00
Concentration technical limits-Upper	700 (675) ¹
SERUM	
Normal volume [μ L]	12
Normal dilution (factor)	1
Below normal volume [μ L]	20
Below normal dilution (factor)	1
Above normal volume [μ L]	2
Above normal dilution (factor)	1
URIN	
Normal volume [μ L]	12
Normal dilution (factor)	1
Below normal volume [μ L]	20
Below normal dilution (factor)	1
Above normal volume [μ L]	2
Above normal dilution (factor)	1
PLASMA	
Normal volume [μ L]	12
Normal dilution (factor)	1
Below normal volume [μ L]	20
Below normal dilution (factor)	1
Above normal volume [μ L]	2
Above normal dilution (factor)	1
CSF	
Normal volume [μ L]	12
Normal dilution (factor)	1
Below normal volume [μ L]	20
Below normal dilution (factor)	1
Above normal volume [μ L]	2
Above normal dilution (factor)	1

¹() With p5p

Results	
Decimals	1
Units	U/L
Correlation factor-Offset	0.000
Correlation factor-Slope	1.000

Range	
Genre	Male
Age	
SERUM	>= <=35.0
URINE	
PLASMA	>= <=35.0
CSF	
Genre	Female
Age	
SERUM	>= <=31.0
URINE	
PLASMA	>= <=31.0
CSF	

Contaminants	
Contaminant 1	
Wash with	
Cycle	
Volume [μ L]	
Contaminant 2	
Wash with	
Cycle	
Volume [μ L]	

Calibrators details	
Calibrator list	Concentration
Cal. 1	0
Cal. 2	*
Cal. 3	*
Cal. 4	*
Cal. 5	*
Cal. 6	*
	Max delta abs.
Cal. 1	0.015
Cal. 2	0.005
Cal. 3	
Cal. 4	
Cal. 5	
Cal. 6	
Drift limit [%]	0.8
Calculations	
Model	X degree
Degree	1

* Enter calibrator value