

# QDx Instacheck™ Anti-CCP Plus

## INTENDED USE

**QDx Instacheck™ Anti-CCP Plus** is a fluorescence Immunoassay (FIA) for the qualitative or semi-quantitative determination of human IgG autoantibodies to cyclic citrullinated peptides (CCP) in human whole blood/serum/plasma. It is useful as an aid in the diagnosis of rheumatoid arthritis (RA) in combination with other clinical and laboratory findings.

For *in vitro* diagnostic use only.

## INTRODUCTION

Rheumatoid Arthritis (RA) is a common, systemic autoimmune disease affecting 0.5-1.0% of the world population. RA is characterized by chronic inflammation of the synovium which can lead to progressive joint destruction, disability and mortality.<sup>(1)</sup> As joint damage is irreversible, early therapeutic intervention is of paramount importance for the prognosis of patients<sup>(2,3)</sup>.

The diagnosis of rheumatological disease are the medical history, clinical findings (including imaging techniques) and serological laboratory tests. Serological diagnostic testing is of growing importance in the early detection and differentiation of RA. The most frequent serological diagnostic testing is the measurement of rheumatoid factor (RF).<sup>(4)</sup> The RF antibody is present in about 75% of RA patients, but its specificity is limited, as it is often present in healthy individuals and patients with other rheumatic or inflammatory diseases, autoimmune diseases or chronic infections.<sup>(5)</sup>

More recently, new specific autoantibodies to citrullinated proteins antigens (ACPAs) have made a crucial contribution to the diagnosis of RA.<sup>(6)</sup> Although many assays are available to test for ACPAs to specific antigens, for the clinical management of RA, most ACPA testing is performed using a synthetic cyclic citrullinated protein (CCP) as the antigen to detect ACPAs. An anti-CCP assay is capable to detect the autoantibodies against citrullinated proteins which have a relatively high sensitivity (reportedly between 50-75%) for rheumatoid arthritis and extremely high specificity (about 90%) for RA.<sup>(7)</sup> Its high specificity is why the anti-CCP test has become an important part of the diagnostic process for RA.

## PRINCIPLE

The test uses a sandwich immunodetection method; Dried antibodies and synthetic cyclic citrullinated peptides (CCPs) in the detector tube, once diluted with the diluent, bind with anti-CCP antibodies in sample to form peptide antigen-antibody complexes. These complexes then migrate through the nitrocellulose matrix and are captured by immobilized streptavidin on the test line. The more anti-CCP antibodies in sample forms the more the peptide antigen/anti-CCP antibodies complex and leads to stronger intensity of fluorescence signal on detector anti-human IgG, which is processed by instrument for **QDx Instacheck** tests to show anti-CCP level in sample.

## COMPONENTS AND REAGENTS

**QDx Instacheck™ Anti-CCP** consists of 'Cartridges', 'Detectors', 'Diluent', 'Capillary tubes' and an 'ID chip'.

- The cartridge contains a test strip, the membrane which has streptavidin at the test line, while chicken IgY at the control line.
- Each cartridge is individually sealed in an aluminum foil pouch containing of a desiccant. 25 sealed cartridges are packed in a box which also contains an ID chip and 25 sealed capillary tubes.
- The detector contains anti-human Immunoglobulin G-

fluorescence conjugate, anti-chicken IgY-fluorescence conjugate, CCP-biotin conjugate, bovine serum albumin (BSA) as a stabilizer and sodium azide in phosphate buffered saline (PBS) as a preservative.

- Each detector contains a granule. 25 tubes of detector are packaged in a pouch.
- The Diluent is dispensed in a vial.
- Detector tube pouch and diluent vial are packaged in a buffer box. The buffer box is packed in a Styrofoam box with icepack for the shipment.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- Follow the instructions and procedures described in this 'Instruction for use'.
- Use only fresh samples and avoid direct sunlight.
- Lot numbers of all the test components (Cartridge, ID chip and detector) must agree each other.
- Do not interchange the test components between different lots or use the test components after the expiration date, either of which might yield misleading test result(s).
- Do not reuse cartridges or detector. A detector should be used for processing of one sample only. A cartridge should be for testing one sample only.
- The cartridge should remain sealed in its original pouch until just before use. Do not use cartridge, if pouch is damaged or has already been opened.
- Frozen sample should be thawed only once. For shipping, samples must be packed in accordance with the local regulations. Sample with severe hemolysis and/or hyperlipidemia must not be used.
- Allow cartridge, detector, diluent and sample to be at room temperature for approximately 30 minutes before use.
- The instrument for QDx Instacheck™ tests may generate slight vibration during use.
- Used detector tubes, diluent vial, pipette tips, capillary tubes and cartridges should be handled carefully and discarded by an appropriate method in accordance with relevant local regulations.
- An exposure to larger quantities of sodium azide may cause certain health issues like convulsions, low blood pressure and heart rate, loss of consciousness, lung injury and respiratory failure.
- **QDx Instacheck™ Anti-CCP Plus** will provide accurate and reliable results subject to the below conditions.

- **QDx Instacheck™ Anti-CCP Plus** should be used only in conjunction with instrument for QDx Instacheck™ tests.

- Have to use recommended anticoagulant sample.

Recommended anticoagulant
Na <sub>2</sub> EDTA, K <sub>2</sub> EDTA, Sodium citrate

## STORAGE AND STABILITY

Component	Storage condition		
	Storage Temperature	Shelf life	Note
Cartridge	4 - 30 °C	20 months	Disposable
Detector tube	2 - 8 °C	20 months	Disposable
Detector diluent	2 - 8 °C	20 months	Unopened
	2 - 8 °C	20 months	opened

- After the cartridge pouch is opened, the test should be performed immediately.

## LIMITATIONS OF THE TEST SYSTEM

- The test may yield false positive result(s) due to the cross-reactions and/or non-specific adhesion of certain sample components to the capture/detector antibodies.
- The test may yield false negative result. The non-responsiveness of the antigen to the antibodies is most common where the epitope is masked by some unknown components, so as not to be detected or captured by the antibodies. The instability or degradation of the antigen with time and/or temperature may cause the false negative as it makes antigen unrecognizable by the antibodies.
- Other factors may interfere with the test and cause erroneous results, such as technical/procedural errors, degradation of the test components/reagents or presence of interfering substances in the test samples.
- Any clinical diagnosis based on the test result must be supported by a comprehensive judgment of the concerned physician including clinical symptoms and other relevant test results.

## MATERIALS SUPPLIED

**REF** IFPC-35

### Components of QDx Instacheck™ Anti-CCP Plus

- Cartridge Box:
  - Cartridges 25
  - 5 µL Capillary tube 25
  - ID Chip 1
  - Instruction for Use 1
- Buffer box:
  - Detector tube 25
  - Detector diluent 1

## MATERIALS REQUIRED BUT SUPPLIED ON DEMAND

Following items can be purchased separately from QDx Instacheck™ Anti-CCP Plus. Please contact our sales division for more information.

- QDx Instacheck™ II **REF** FPRR021
- Boditech Anti-CCP Plus Control **REF** CFPO-288
- i-chamber **REF** FPRR009

## SAMPLE COLLECTION AND PROCESSING

The sample type for QDx Instacheck™ Anti-CCP Plus is human whole blood/serum/plasma.

- It is recommended to test the sample within 24 hours after collection.
- The serum or plasma should be separated from the clot by centrifugation within 3 hours after the collection of whole blood.
- The serum and plasma, whole blood may be stored for up to a month at 2-8 °C prior to being tested. If testing will be delayed more than a month, samples should be frozen at -20 °C.
- The freezing storage of sample up to 3 months does not affect the quality of results. However, the whole blood sample should not be kept in a freezer in any case.
- Once the sample was frozen, it should be used one time only for test, because repeated freezing and thawing can result in the change of test values.
- Fingertip blood sample should be collected as follows:

- Position the hand with the palm facing upwards. Blood should be normally drawn from the middle or ring finger of the non-dominant hand. Apply intermittent pressure towards its tip.
- Wipe the fingertip clean with an alcohol pad.
- Allow the finger to dry completely because blood will not form a drop if the puncture site is moist and because the residual alcohol at the fingertip may dilute the blood sample and affect the test result.
- Hold the finger and puncture the fingertip by firmly pressing a new sterile lancet against it at an off-center position.
- Wipe away the first drop of blood with a sterile gauze pad or cotton ball.
- Massage the finger towards its tip to form a new drop of blood. Blood will flow easily if the finger is held lower than the elbow.
- Hold the handle of a capillary tube and touch the mouth of the capillary to the drop of blood.
- Let the blood fill the capillary tube completely.
- It may be sometimes necessary to massage the finger again for an additional drop of blood for filling the capillary tube.

## TEST SETUP

- Check the components of QDx Instacheck™ Anti-CCP Plus: Sealed Cartridge, Detector Tube, Diluent vial, capillary tube and ID Chip.
- Ensure that the lot number of the cartridge matches with that of the ID chip as well as the buffer box.
- Keep the sealed cartridge (if stored in refrigerator) and the detection buffer tube at room temperature for at least 30 minutes just prior to performing the test. Place the cartridge on a clean, dust-free and flat surface.
- Turn on the instrument QDx Instacheck™ tests. (Please refer to the 'Instrument for Operation Manual' for complete information and operating instructions.)

## Caution

- To minimize erroneous test results, we suggest that the ambient temperature of the cartridge should be 25 °C during the reaction time after loading sample mixture to the cartridge
- To maintain the ambient temperature to 25 °C, you can use various devices such as an i-Chamber or an incubator and so on.

## TEST PROCEDURE

1. Transfer 150 µL of diluent using a pipette to a tube containing detector.
2. Transfer 5 µL of sample (Human whole blood/serum/plasma/control) to the detector tube.
  - \* If you use a capillary tube (5 µL), put it into the detector tube after collecting whole blood sample.
3. Close the lid of the detector tube and mix the sample thoroughly by shaking it about 20 times.
4. Pipette out 75 µL of a sample mixture and load it into the sample well on the cartridge.
5. Insert and leave the sample mixture-loaded cartridge in the i-Chamber or incubator (25 °C) for 12 minutes.
6. To scan the sample mixture-loaded cartridge, insert it into the cartridge holder of the instrument for QDx Instacheck™ tests. Ensure proper orientation of the cartridge before pushing it all the way inside the cartridge holder. An arrow has been marked on the cartridge especially for this purpose.
  - ⚠ Scan the sample mixture-loaded cartridge immediately when the incubation time is over. If not, it will cause inexact test result.**

7. Tap the 'Start' button on the instrument for QDx Instacheck™ tests to start the scanning process
  8. The instrument for QDx Instacheck™ tests will start scanning the sample-loaded cartridge immediately.
  9. Read the test result on the display screen of the QDx Instacheck™ tests.
- ※ Please refer to the 'Instrument for QDx Instacheck™ tests Operation Manual' for complete information and operating instructions.

### INTERPRETATION OF TEST RESULT

- The instrument for QDx Instacheck™ tests calculates the test result automatically and displays anti-CCP concentration and anti-CCP state of the test sample.

Test result [U/mL]	Display [U/mL]
< 5.0	< 3.5 or value, (Neg)
5.0 ≤, < 300	Value, (Pos)

- The cut-off: 5 U/mL
- Working range: 3.5-300 U/mL

### QUALITY CONTROL

- Quality control tests are a part of the good testing practice to confirm the expected results and validity of the assay and should be performed at regular intervals.
- The control tests should be performed immediately after opening a new test lot to ensure the test performance is not altered.
- Quality control tests should also be performed whenever there is any question concerning the validity of the test results.
- Control materials are not provided with **QDx Instacheck™ Anti-CCP Plus**. For more information regarding obtaining the control materials, contact the technical section at Diasys Diagnostics India Private Limited.  
(Please refer to the instruction for use of control material.)

### PERFORMANCE CHARACTERISTICS

#### Analytical sensitivity

Limit of Blank (LoB):	2.32 IU/mL
Limit of Detection (LoD):	3.49 IU/mL

#### Analytical Specificity

##### Interference

There was no significant interference from these materials with the **QDx Instacheck™ Anti-CCP Plus** test measurements.

Interference material	Concentration
Hemoglobin	500 mg/dL
Bilirubin	0.2 mg/mL
Triglyceride	2000 mg/dL
Rheumatoid factor	78 IU/mL
Human serum albumin	12 g/dL

##### Cross-reactivity test

There was no significant cross-reactivity from these materials with the **QDx Instacheck™ Anti-CCP Plus** test measurements. (α-SSA, α-SSB, α-Sm, α-RNP, α-ds-DNA, α-Jo-1, α-Scl-70, α-Ribo-P, anti-nuclear antibody(ANA))

#### Precision

##### Between lot

One person tested three different lots of **QDx Instacheck™ Anti-CCP Plus**, ten times at each concentration of the control standard.

##### Between person

Three different persons tested **QDx Instacheck™ Anti-CCP Plus**, ten times at each concentration of the control standard.

##### Between day

One person tested **QDx Instacheck™ Anti-CCP Plus** for five days, five times at each concentration of the control standard.

##### Between site

One person tested **QDx Instacheck™ Anti-CCP Plus** at three different sites, five times at each concentration of the control standard.

Anti-CCP [U/mL]	Between lot		Between person		Between day		Between site	
	AVG	CV(%)	AVG	CV(%)	AVG	CV(%)	AVG	CV(%)
6.25	6.22	5.0	6.21	5.9	6.25	6.6	6.13	5.2
30.00	29.61	5.7	30.39	5.3	29.98	5.7	29.50	6.2
100.00	99.42	6.6	99.27	6.5	100.64	5.6	98.55	6.9

#### Accuracy

The accuracy was determined by 3 different lots testing six times each human serum.

Anti-CCP [U/mL]	Lot 1	Lot 2	Lot 3	Mean	CV (%)	Bias (%)
4.15	3.83	3.86	3.87	3.86	2.5	-7.1
5.31	5.05	4.88	4.91	4.95	3.3	-6.9
16.66	15.21	15.32	15.17	15.23	1.8	-8.6
19.79	18.53	18.07	18.38	18.33	4.0	-7.4
67.02	67.16	67.67	66.55	67.13	2.2	0.2

#### Comparability













Total (N=216)	n	QDx Instacheck™ Anti-CCP Plus	
		Positive	Negative
<b>Axis-shield</b>	Pos. 116	109	7
<b>FCCP600</b>	Neg. 100	4	96
Total sample		216	
Positive Agreement rate (≥5 U/mL) (%)		93.9	
Negative Agreement rate (≥5 U/mL) (%)		96.0	
Total (%)		94.9	

### REFERENCES

1. Birch JT Jr, Bhattacharya S. Emerging trends in diagnosis and treatment of rheumatoid arthritis. Prim Care. 2010;37:779–92.
2. Landewe RB. The benefits of early treatment in rheumatoid arthritis: confounding by indication, and the issue of timing. Arthritis Rheum 2003;48(1):1-5.
3. Lard LR, Visser H, Speyer I, et al. Early versus delayed treatment in patients with recent-onset rheumatoid arthritis: comparison of two cohorts who received different treatment strategies. Am J Med 2001;111:446-51.
4. Mageed RA. The RF antigen. In Manual of Biological Markers of Disease Edited by van Venrooij WJ & Maini RN Dordrecht: Kluwer Academic Publishers; 1996. pp. 1–27.
5. van Boekel MA, Vossenaar ER, van den Hoogen FH, van Venrooij WJ. Autoantibody systems in rheumatoid arthritis: specificity, sensitivity and diagnostic value. Arthritis Res 2002;4:87–93
6. Systematic review: accuracy of anti-citrullinated Peptide antibodies for diagnosing rheumatoid arthritis. Whiting PF, Smidt N, Sterne JA, Harbord R, Burton A, Burke M, Beynon R, Ben-Shlomo Y, Axford J, Dieppe P Ann Intern Med. 2010 Apr 6; 152(7):456-64; W155-66.
7. The diagnostic properties of rheumatoid arthritis antibodies

recognizing a cyclic citrullinated peptide. Schellekens GA, Visser H, de Jong BA, van den Hoogen FH, Hazes JM, Breedveld FC, van Venrooij WJ Arthritis Rheum. 2000 Jan; 43(1):155-63.

**Note:** Please refer to the table below to identify various symbols.

	Sufficient for <n> tests
	Read instruction for use
	Use by Date
	Batch code
	Catalog number
	Caution
	Manufacturer
	Authorized representative of the European Community
	In vitro diagnostic medical device
	Temperature limit
	Do not reuse
	This product fulfills the requirements of the Directive 98/79/EC on in vitro diagnostic medical devices



#### **Boditech Med Incorporated**

43, Geodudanji 1-gil, Dongnae-myeon,  
Chuncheon-si, Gang-won-do, 24398  
Republic of Korea  
Tel: +(82) -33-243-1400  
Fax: +(82) -33-243-9373  
www.boditech.co.kr



#### **Obelis s.a**

Bd. Général Wahis 53,  
1030 Brussels, BELGIUM  
Tel: +(32) -2-732-59-54  
Fax: +(32) -2-732-60-03  
E-Mail: mail@obelis.net

Imported and Marketed by  
**DiaSys Diagnostics India Private Limited**  
Plot no. A-821, TTC Industrial Area,  
MIDC, Mahape,  
Navi Mumbai, 400710, Maharashtra, India

For feedback/queries contact customer care at:  
Toll Free number: 18001201447

Email ID: [Helpdesk.Service@diasys.in](mailto:Helpdesk.Service@diasys.in)  
Website: [www.diasys.in](http://www.diasys.in)

Revision No. 00

Date of last revision: January 22, 2021

